





REORDER NUMBER:

RTG-01006

RespVent™ HME (HCH) Non-Filtered Hygroscopic Condenser Humidifier With Gas Sampling Port



Clean, Ready to Use
For Single Patient Use

RxOnly

Caution: Federal law (USA) restricts this device to use by or at the direction of a physician

INSTRUCTIONS FOR USE

Read these instructions carefully before using the product. Place the *RespVent™ HMEF(HCH)* between the proximal end of the artificial airway and the Y-piece of the breathing circuit. Always replace the *RespVent™ HMEF(HCH)* after each patient. When used continuously on a single patient, change the *RespVent™ HMEF(HCH)* every 24 hours or more frequently as required.

Warning: The *RespVent™ HMEF(HCH)* is designed for single patient use only and must not be cleaned and reused.

CONTRAINDICATIONS

The *RespVent™ HMEF(HCH)* is contraindicated in patients producing fulminating, frothy secretions within their airway and lungs. The *RespVent™ HMEF(HCH)* shall not be used on patients with very small tidal volumes, for example, neonates. The *RespVent™ HMEF(HCH)* shall not be used together with active humidifier or nebulizers.

PRECAUTIONS

All tubing and connections to the *RespVent™ HMEF(HCH)* shall be properly attached and checked for leakage prior to use. Compensation of ventilation may be necessary when using the *RespVent™ HMEF(HCH)* since dead space will be added to the system. During the use of the *RespVent™ HMEF(HCH)*, the patient shall be closely monitored and proper airway care administered if complications arise. The *RespVent™ HMEF(HCH)* must be changed between patients.

| | | |
|--------------------------|--|----------------------------|
| Recommended tidal volume | 150-1500ml | |
| Pressure drop | at flow 30L/min 1.0 cm H ₂ O at flow 60L/min 2.8 cm H ₂ O | |
| Internal volume | 34ml | |
| Moisture loss | (per ISO 9360 test procedure) | |
| | Tidal Volume = 250 ml | 6.1 mg H ₂ O/L |
| | Tidal Volume = 500 ml | 8.3 mg H ₂ O/L |
| | Tidal Volume = 1000 ml | 10.3 mg H ₂ O/L |

Weight 15 g
Connections 22M/15F - 15M

STORAGE TEMPERATURE:
-30°C/-22°F
+40°C/+104°F

For more info

Telephone
323-761-0035

E-mail
info@rtghealth.com

www.rtghealth.com