

REORDER NUMBER:


# RTG-01020 RespVent<sup>™</sup> Filtered (HCH)

Small-Volume Filtered Hygroscopic Condenser **Humidifier With Gas Sampling Port** 

## **Pediatric**

## RxOnly

Caution: Federal law (USA) restricts this device to use by or at the direction of a physician













Clean. Ready to Use For Single Patient Use

### **INSTRUCTIONS FOR USE**

Read these instructions carefully before using the product. Place the *RespVent™ HMEF(HCH)* between the proximal end of the artificial airway and the Y-piece of the breathing circuit. Always replace the  $RespVent^{m}$  HMEF(HCH) after each patient. When used continuously on a single patient, change the  $RespVent^{TM}$  HMEF(HCH) every 24 hours or more frequently as required. **Warning:** The  $RespVent^{TM}$  HMEF(HCH) is designed for single patient use only and must not be cleaned and reused.

#### CONTRAINDICATIONS

The RespVent™ HMEF(HCH) is contraindicated in patients producing fulminating, frothy secretions within their airway and lungs. The *RespVent<sup>TM</sup> HMEF(HCH)* shall not be used on patients with very small tidal volumes, for example, neonates. The  $RespVent^{TM}$  HMEF(HCH) shall not be used together with active humidifier or nebulizers.

#### **PRECAUTIONS**

All tubing and connections to the  $RespVent^{TM}$  HMEF(HCH) shall be properly attached and checked for leakage prior to use. Compensation of ventilation may be necessary when using the *RespVent™HMEF(HCH)* since dead space will be added to the system. During the use of the *RespVent<sup>TM</sup> HMEF(HCH)*, the patient shall be closely monitored and proper airway care administered if complications arise. The *RespVent™ HMEF(HCH)* must be changed between patients.

Recommended tidal volume Pressure drop

Internal volume Moisture loss

For more info

70~600ml

at flow 15L/min 0.6 cm H<sub>2</sub>O at flow 30L/min 1.6 cm H<sub>2</sub>0

(per ISO 9360 test procedure) Tidal Volume = 250 ml

Tidal Volume = 500 ml

10.5 mg H<sub>2</sub>O/L 16.0 mg H<sub>2</sub>0/L Filtration Efficiency Bacterial(%) Viral(%) Weight Connections

99.999% 99.98% 18 g 22M/15F-15M

STORAGE TEMPERATURE: -30°C/-22°F +40°C/+104°F