

REORDER NUMBER:


RTG-01021 RespVent[™] HME (HCH)

Small-Volume Non-Filtered Hygroscopic Condenser Humidifier With Gas Sampling Port

Pediatric

RxOnly

Caution: Federal law (USA) restricts this device to use by or at the direction of a physician



Clean, Ready to Use For Single Patient Use

INSTRUCTIONS FOR USE

Read these instructions carefully before using the product. Place the $RespVent^{TM}$ HMEF(HCH) between the proximal end of the artificial airway and the Y-piece of the breathing circuit. Always replace the $RespVent^{TM}$ HMEF(HCH) after each patient. When used continuously on a single patient, change the $RespVent^{TM}$ HMEF(HCH) every 24 hours or more frequently as required. **Warning:** The $RespVent^{TM}$ HMEF(HCH) is designed for single patient use only and must not be cleaned and reused.

CONTRAINDICATIONS

The *RespVentTM HMEF(HCH)* is contraindicated in patients producing fulminating, frothy secretions within their airway and lungs. The *RespVentTM HMEF(HCH)* shall not be used on patients with very small tidal volumes, for example, neonates. The *RespVentTM HMEF(HCH)* shall not be used together with active humidifier or nebulizers.

PRECAUTIONS

All tubing and connections to the *RespVentTM HMEF(HCH)* shall be properly attached and checked for leakage prior to use. Compensation of ventilation may be necessary when using the *RespVentTM HMEF(HCH)* since dead space will be added to the system. During the use of the *RespVentTM HMEF(HCH)*, the patient shall be closely monitored and proper airway care administered if complications arise. The *RespVentTM HMEF(HCH)* must be changed between patients.

Recommended tidal volume Pressure drop

Internal volume Moisture loss 80~600ml

at flow 30L/min $\,$ 0.8 cm $\,$ H20 at flow 60L/min $\,$ 2.2 cm H20 $\,$

21mi

(per ISO 9360 test procedure)

Tidal Volume = 250 ml Tidal Volume = 500 ml Weight Connections 9 g

22M/15F-15M

STORAGE TEMPERATURE: -30°C/-22°F +40°C/+104°F



9.3 mg H₂0/L

11.3 mg H₂0/L